# ImmuLisa™ Anti-Cardiolipin Antibody (ACA) Screen ELISA

For *in vitro* Diagnostic Use

CLIA Complexity: High  
CDC Analyte Identification Code: 0434  
CDC Test System Identification Code: 28515

## PRODUCT INSERT

Catalog No. 1118S

**ACA Screen**

96 Determinations

## INTENDED USE

An enzyme linked immunoassay (ELISA) for the detection of IgG, IgA and IgM class anticardiolipin antibodies as an aid in assessing the risk of thrombosis in individuals with anti-phospholipid syndrome.

## SUMMARY AND EXPLANATION

Antiphospholipid antibodies are a heterogeneous group of autoantibodies against negatively charged phospholipids<sup>1</sup>. They are detected primarily by the anti-cardiolipin antibody (ACA) test, the biological false positive test for syphilis and the lupus anticoagulant test. These three tests detect related, but not necessarily identical antibodies. Thus, more than one of these tests is sometimes necessary to identify antiphospholipid antibodies.

Of the various tests for the detection of antiphospholipid antibodies, the anti-cardiolipin antibody test performed by ELISA is the most sensitive<sup>2</sup>. The presence of anti-cardiolipin antibodies helps to identify patients at risk of venous and/or arterial thrombosis often accompanied by thrombocytopenia, a syndrome referred to as antiphospholipid syndrome<sup>1-12</sup>. The antiphospholipid syndrome most commonly occurs in patients with systemic lupus erythematosus (SLE) or lupus-like disease where the criteria for SLE are not fulfilled<sup>5-7</sup>. High levels of anti-cardiolipin antibodies occur in thrombosis, fetal loss, thrombocytopenia and several other disorders<sup>1-15</sup>. Low levels of anti-cardiolipin antibodies are found in a variety of clinical disorders which are unrelated to antiphospholipid syndrome. Therefore, low levels of these antibodies are of limited significance.

IgG and IgA class anti-cardiolipin antibodies appear to be more closely associated with antiphospholipid syndrome than the IgM class antibodies. However, IgM antibodies appear to be more influenced by treatment<sup>5,10</sup>. Low levels of IgM antibodies can be identified in other autoimmune diseases such as rheumatoid arthritis, primary Sjögren's Syndrome, drug induced lupus erythematosus, Lyme disease, and syphilis<sup>8,10</sup>.

The ImmuLisa™ ACA Screen ELISA offers a convenient method of screening patient sera for the presence of ACA before further identifying positive samples for individual isotypes with the ImmuLisa™ ACA-IgG, -IgA and -IgM ELISA.

## PRINCIPALS OF THE PROCEDURE

Cardiolipin antigen is bound to the wells of a polystyrene microwell plate followed by blocking the unreacted sites to reduce non-specific binding. Controls, calibrator and diluted patient sera are added to separate wells, allowing any cardiolipin antibodies present to bind to the immobilized antigen. Unbound sample is washed away and an enzyme labeled anti-human conjugate is added to each well. These enzyme conjugated antibodies bind specifi-

cally to the human immunoglobulin of the appropriate class. After washing away any unbound conjugate, specific enzyme substrate (pNPP) is then added to the wells. After stopping the enzymatic reaction, the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in ELISA units per milliliter (EU/ml).

**REAGENTS**

**Storage and Preparation**

Store all reagents at 2-8°C. **Do not freeze.** Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water. Coated microwell strips are for one time use only.

**Precautions**

For *in vitro* Diagnostic Use. All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>17</sup>. **WARNING** - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources other than the same catalog number from IMMCO DIAGNOSTICS.

Use good laboratory techniques to minimize microbial and chemical contamination. Do not use after expiration date.

**Materials provided**

ImmuLisa™ ACA Screen ELISA REF 1118S

Kits contain sufficient reagents to perform 96 determinations each.

12 x 8	MICROPLATE ACA	<b>Microplate</b> with individual breakaway microwells coated with cardiolipin antigen.
1 x 1.5 ml	CALIBRATOR ACA *	Ready to use <b>Calibrator A</b> ( <i>green cap</i> ). Human serum containing antibodies to cardiolipin.
1 x 1.5 ml	CONTROL + ACA *	Ready to use <b>Positive Control</b> ( <i>red cap</i> ). Contains human serum positive for ACA.
1 x 1.5 ml	CONTROL - *	Ready to use <b>Negative Control</b> ( <i>white cap</i> ). Contains human serum.
1 x 12 ml	CONJ ALKPHOS *	Ready to use <b>anti-human Alk. Phos. Conjugate</b> . Color coded pink.
2 x 60 ml	DIL *	Ready to use <b>Serum Diluent</b> . Color coded blue.

1 x 12 ml **SUBSTRATE**\*

Ready to use **Enzyme Substrate**.  
Contains pNPP. **Protect from light.**

1 x 12 ml **STOP**

Ready to use **Stop Solution**.

2 **BUF WASH**

Powder **Wash Buffer**. Reconstitute  
to one liter each.

\* Contains < 0.1% NaN<sub>3</sub>

### Materials Required But Not Provided

- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Deionized or distilled water
- Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.
- Squeeze bottle to hold diluted wash buffer
- Timer
- Absorbent paper
- Automatic microplate washer capable of dispensing 200 µl

### SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

### PROCEDURE

#### Procedural Notes

- Before starting the assay read these instructions carefully.
- Bring all reagents and samples to room temperature (20-26°C) for 30 minutes. Return materials to refrigerator immediately after use.
- Prepare all dilutions of the patient specimens before starting the test.
- **Immediately return unused strips to the pouch containing desiccants and seal securely to minimize exposure to water vapor.**
- Wash step: Good technique is critical. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- Careful timing is important. Incubation periods begin after dispensing reagents.

#### Assay procedure

1. **ALL REAGENTS MUST BE BROUGHT TO ROOM TEMPERATURE (20-26°C) PRIOR TO BEGINNING THE ASSAY.**
2. Label protocol record to indicate specimen placement in the microplate. It is good laboratory practice to test specimens in duplicate.
3. Prepare a **1:201** dilution of the patient specimen by mixing **5 µl** of the patient specimen with **1.0 ml** of Serum Diluent.
4. Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.

## QUALITATIVE DETERMINATION

A	BLANK	S5		
B	NEG	S6		
C	POS	S7		
D	CAL D	S8		
E	S1	S9		
F	S2	S10		
G	S3	S11		
H	S4	S12		
	1	2	3	4

5. Add **100 µl** of Calibrators, Positive and Negative controls and diluted patient specimens to the appropriate microwells indicated on the protocol record.
- Note:** Include one well with **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank. The absorbance of this well should not be greater than 0.3.
6. Incubate **30 minutes** ( $\pm$  5 minutes) at room temperature on a level surface.
7. Wash step: Thoroughly aspirate the contents of each well. Add 200-300µL of the **reconstituted** wash buffer to all wells then aspirate. Repeat this sequence thrice more for a total of four washes. Invert the plate and tap it on absorbent material to remove any residual fluid after the last wash. Do not dry wells completely.
8. Add 100µL of the Conjugate to each well.
9. Incubate the wells for 30 minutes at room temperature.
10. Wash step: Repeat step 7.
11. Add 100µL of Enzyme Substrate to each well.
12. Incubate for 30 minutes at room temperature.
13. Add 100µL of Stop Solution to each well. Maintain the same sequence and timing of Stop Solution addition as was used for the Enzyme Substrate. Read the absorbance (OD) of each well at 405nm within one hour of stopping the reaction.
14. Read the absorbance (OD) of each well at 405nm using a single or dual wavelength microplate reader against the reagent blank set at zero absorbance.

### Quality Control

Calibrator, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be  $<0.3$ . The Calibrator should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be  $<20$  EU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining EU/ml.

### RESULTS

#### Calculations

The following method should be used to determine whether the specimen is positive or negative for ACA.

## QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

### Interpretation

The following information serves only as a guide in the interpretation of the laboratory results. Each laboratory must determine its own normal values.

ACA Screen Value Value/Valeur/Valore/Valor/Wert	Interpretation/Interprétation/Interpretazione/ Interpretación/Deutung/Interpretação
<20 EU/ml	Neg (-)
20-25 EU/ml	Borderline/Indéterminé/Indeterminato/Indeterminado/ Unbestimmt
>25 EU/ml	Pos (+)

The literature suggests that, low positive anti-cardiolipin antibody levels may occur in a variety of clinical disorders unrelated to antiphospholipid antibody syndrome. Hence according to the investigators recommendations the diagnosis of antiphospholipid antibody syndrome should be made only when the test results are moderately or highly positive<sup>14</sup>.

### LIMITATIONS OF THE PROCEDURE

A diagnosis cannot be made on the basis of ACA results alone. The results of other laboratory tests and clinical findings must be considered. Tests for rheumatoid factor (RF) should be performed because RF may interfere with this assay. When a negative ACA test occurs in the presence of clinical indications, a lupus anticoagulant test or other additional testing is indicated. ACA also occur transiently in a variety of infectious diseases. In these cases patients positive for ACA should be retested following an appropriate interval. Confirmed active or seropositive syphilis patients can have elevated ACA levels. To rule out syphilis, confirmatory tests should be performed.

### EXPECTED VALUES

The incidence of ACA in various disease conditions is summarized in the following tables:

#### *Incidence of ACA in SLE*<sup>15,18</sup>

Antibody Isotype	% Incidence
IgG	39-44
IgA	17-57
IgM	5-33
Any isotype	53

#### *Disease Association of ACA*

Condition	% Incidence
Recurrent Venous Thrombosis	28-71
Recurrent Fetal Loss	28-64
Transverse Myelitis	50
Hemolytic Anemia	38
Thrombocytopenia	27-33
Arterial Occlusions	25-31
Livedo Reticularis	25
Pulmonary Hypertension	20-40

## PERFORMANCE CHARACTERISTICS

### Comparison of ELISA Methods for the Detection of ACA

		IMMCO IgG-ACA		
		POS (+)	NEG (-)	TOT (=)
ImmuLisa™	POS (+)	31	0	31
ACA IgG/A/M	NEG (-)	2	42	44
	TOT (=)	33	42	75

relative specificity/spécificité relative/especificidad relativa/relative Spezifität/specificità relativa/especificidade relativa: 95.5%

relative sensitivity/sensibilité relative/sensibilidad relativa/relative Sensitivität/sensibilità relativa/sensibilidade relativa: 100.0%

relative agreement/concordance/concordancia relativa/relative Übereinstimmung/concordanza relativa/acordo relativo: 97.3%

### Precision:

Two samples with known concentrations of ACA were assayed in 10 replicates over a period of two weeks. Intra-and inter-assay coefficient of variation (CV) were as follows:

	Intra- plate plaque/piastra/placa/Platte CV	Inter-plate plaque/piastra/placa/Platte CV
1.	7.6%	5.9%
2.	10.4%	6.3%

### Recovery:

Samples with known concentrations of ACA were mixed with the appropriate dilutions of another positive sample with known amounts of ACA. The ACA values of the spiked samples were then determined and from the values obtained the percent recovery calculated.

	ACA added/supplémentaire/ aggiunto/agregado/ hinzugefügt/adicionado	ACA obtained/obtenu/ verificato/obtenido/ erreicht/obtido	% Recovery/Rétablissement/ Recupero/Recuperación/ Wiederaufnahme/Recuperação
1.	127.0	135.8	106.9
2.	103.8	103.2	99.4
3.	64.3	61.2	95.2

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<b>Doença</b>	<b>% Incidência</b>
Recurrent Venous Thrombosis	28-71
Recurrent Fetal Loss	28-64
Transverse Myelitis	50
Hemolytic Anemia	38
Thrombocytopenia	27-33
Arterial Occlusions	25-31
Livedo Reticularis	25
Pulmonary Hypertension	20-40

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## REFERENCES•REFERENCIAS•LITERATUR•RIFERIMENTI

1. Harris EN, Phil M, Asherson RA et al. Antiphospholipid antibodies - antibodies with a difference. *Ann Rev Med*; 1988, 39:261-271.
2. Harris EN, Gharavi AE, Wasley GD and Hughes GRV. Use of an enzyme-linked immunosorbent assay and inhibition studies to distinguish between antibodies to cardiolipin from patients with syphilis or autoimmune disorders. *J Infect Dis*; 1988, 157:23-31.
3. Lockshin MD. Anti-cardiolipin antibody. *Arth Rheum*; 1987, 30:471-472.
4. Derksen RHWM, Beisma D, Bouma BN et al. Discordant effects of prednisone on anticardiolipin antibodies and the lupus anti-coagulant. *Arth Rheum*; 1986, 29:1295-1296.
5. Alarcón-Segovia D, Deleze M, Oria CV et al. Antiphospholipid antibodies and the antiphospholipid syndrome in systemic lupus erythematosus. *Medicine*; 1989, 68:353-365.
6. Triplett DA. Antiphospholipid antibodies and recurrent pregnancy loss. *Am J Reprod Immunol*; 1989, 20:52-67.
7. Montecucco C, Longhi M, Caporali R et al. Hematological abnormalities associated with anticardiolipin antibodies. *Hematologica*; 1989, 74:195-204.
8. Asherson RA, Khamashta MA, ORDI-Ros J et al. The "primary" antiphospholipid syndrome: major clinical and serological features. *Medicine*; 1989, 68:366-374.
9. Walport MJ. Pregnancy and antibodies to phospholipids. *Ann Rheum Dis*; 1989, 48:795-797.
10. Machworth-Young CG, Loizou S and Walport MJ. Antiphospholipid antibodies and disease. *Medicine*; 1989, 72:767-777.
11. Harris EN, Gharavi AE and Hughes GRV. Antiphospholipid antibodies. *Clin Rheumatol*; 1985, 1:591-609.
12. Lubbe WF, Palmer SJ, Butler WS et al. Fetal survival after prednisone suppression of maternal lupus anticoagulant. *Lancet*; 1983, 1:1361-1363.
13. Hogan MJ, Brunet DG, Ford PM et al. Lupus anticoagulant, antiphospholipid antibodies and migraine. *Can J Neurological Sci*; 1988, 15:420-425.
14. Harris EN and Pierangeli S. Antiphospholipid antibodies: method of detection. *Am J Reproduct Immunol*; 1992, 28:208-210.
15. Silveira LH, Lopez LR, Uzzle RT et al. IgA, IgG and IgM anticardiolipin antibodies in black patients with systemic lupus erythematosus. *Arthr Rheum 56th Annual Meeting*; 1992, S125.
16. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health; 1993, (HHS Pub. No [CDC] 93-8395).
17. Alarcón-Segovia D, Delezé M, Oria CV et al. Antiphospholipid antibodies and the antiphospholipid syndrome in systemic lupus erythematosus. A perspective analysis of 500 consecutive patients. *Medicine*; 1987, 68:353-365.

**NOTES:**

**NOTES:**

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**

**60 Pineview Drive**

**Buffalo, NY 14228-2120**

**Telephone: (716) 691-0091**

**Fax: (716) 691-0466**

**Toll Free USA/Canada: 1-800-537-TEST**

**E-Mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,

The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

[www.emergogroup.com](http://www.emergogroup.com)