# Immuli<sup>TM</sup> Rheumatoid Factor IgM ELISA

IVD

## PRODUCT INSERT

REF	Code	1139	96 Determinations
-----	------	------	-------------------

## INTENDED USE

An enzyme linked immunosorbant assay (ELISA) for the detection and quantitation of Rheumatoid Factor (RF) IgM isotype in human serum.

## SUMMARY AND EXPLANATION

The measurement of RF is important in the diagnosis and prognosis of rheumatoid arthritis, as high titers of RF occur in sera of patients who tend to develop extra-articular complications<sup>1,2</sup>. The majority of routine laboratory tests measure IgM RF by its ability to agglutinate sheep red blood cells, latex or similar particles coated with IgG<sup>1-4</sup>.

RF is present in 70-90% of patients with rheumatoid arthritis and is included in the classification criteria<sup>5</sup>. According to the revised ARA criteria, if RF is positive in patients with arthritis of three or more joints then the patient has rheumatoid arthritis. Arthritis of less than three joints with RF negative excludes rheumatoid arthritis. Such a classification affords 93.5% sensitivity and 89.3% specificity for rheumatoid arthritis<sup>5</sup>.

RF as detected by agglutination is of the IgM isotype. Other methods such as ELISA have improved specificity, sensitivity and reliability over existing and routinely used agglutination methods<sup>3,4,6</sup>. ELISA methods can detect RF of various immunoglobulin isotypes. Such a distinction is not possible with traditional agglutination tests.

## PRINCIPALS OF THE PROCEDURE

IgG is bound to the wells of a polystyrene microwell plate followed by blocking the unreacted sites to reduce non-specific binding. Controls, calibrators and diluted patient sera are added to separate wells, allowing any RF present to bind. Unbound sample is washed away and an enzyme labeled anti-human conjugate is added to each well. These enzyme conjugated antibodies bind specifically to the human immunoglobulin. After washing away any unbound conjugate, specific enzyme substrate (pNPP) is then added to the wells. After stopping the enzymatic reaction, the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in International units per milliliter (IU/ml, WHO).

## REAGENTS

### Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water. Coated microwell strips are for one time use only.

### Precautions

For *in vitro* Diagnostic Use. All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>17</sup>. **WARNING** - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly

explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources other than the same catalog number from IMMCO DIAGNOSTICS.

Use good laboratory techniques to minimize microbial and chemical contamination. Do not use after expiration date.

**Materials provided**

ImmuLisa™ IgM-RF ELISA **REF** 1139

Kits contain sufficient reagents to perform 96 determinations each.

<b>12 x 8</b>	<b>MICROPLATE RF</b>	<b>Microplate</b> with individual breakaway microwells coated with IgG antigen.
<b>1 x 1.5 ml</b>	<b>CALIBRATOR A IgM-RF</b> *	Ready to use <b>Calibrator A</b> ( <i>green cap</i> ). Human serum containing RF antibodies.
<b>1 x 1.5 ml</b>	<b>CALIBRATOR B IgM-RF</b> *	Ready to use <b>Calibrator B</b> ( <i>violet cap</i> ). Human serum containing RF antibodies.
<b>1 x 1.5 ml</b>	<b>CALIBRATOR C IgM-RF</b> *	Ready to use <b>Calibrator C</b> ( <i>blue cap</i> ). Human serum containing RF antibodies.
<b>1 x 1.5 ml</b>	<b>CALIBRATOR D IgM-RF</b> *	Ready to use <b>Calibrator D</b> ( <i>yellow cap</i> ). Human serum containing RF antibodies.
<b>1 x 1.5 ml</b>	<b>CONTROL + IgM-RF</b> *	Ready to use <b>Positive Control</b> ( <i>red cap</i> ). Contains human serum positive for RF antibodies.
<b>1 x 1.5 ml</b>	<b>CONTROL -</b> *	Ready to use <b>Negative Control</b> ( <i>white cap</i> ). Contains human serum.
<b>1 x 12 ml</b>	<b>IgM-CONJ ALKPHOS</b> *	Ready to use <b>anti-human Alk. Phos. Conjugate</b> . Color coded pink.
<b>1 x 60 ml</b>	<b>DIL</b> *	Ready to use <b>Serum Diluent</b> . Color coded blue.
<b>1 x 12 ml</b>	<b>SUBSTRATE</b> *	Ready to use <b>Enzyme Substrate</b> . Contains pNPP. <b>Protect from light.</b>
<b>1 x 12 ml</b>	<b>STOP</b>	Ready to use <b>Stop Solution</b> .
<b>2</b>	<b>BUF WASH</b>	Powder <b>Wash Buffer</b> . Reconstitute to one liter each.

\* Contains < 0.1% NaN<sub>3</sub>

**Materials Required But Not Provided**

- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Deionized or distilled water
- Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.
- Squeeze bottle to hold diluted wash buffer

- Timer
- Absorbent paper
- Automatic microplate washer capable of dispensing 200 µl

**SPECIMEN COLLECTION AND HANDLING**

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

**PROCEDURE**



**Procedural Notes**

- Before starting the assay read these instructions carefully.
- Bring all reagents and samples to room temperature (20-26°C) for 30 minutes. Return materials to refrigerator immediately after use.
- Prepare all dilutions of the patient specimens before starting the test.
- **Immediately return unused strips to the pouch containing desiccants and seal securely to minimize exposure to water vapor.**
- Wash step: Good technique is critical. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- Careful timing is important. Incubation periods begin after dispensing reagents.

**Assay procedure**

1. **ALL REAGENTS MUST BE BROUGHT TO ROOM TEMPERATURE (20-26°C) PRIOR TO BEGINNING THE ASSAY.**
2. Label protocol record to indicate specimen placement in the microplate. It is good laboratory practice to test specimens in duplicate.
3. **Qualitative determination:** use only Calibrator D.  
**Semi-quantitative determination:** use Calibrators A - D as shown in the example below.

**Qualitative determination:      Semi-quantitative determination:**

4. Prepare a **1:101** dilution of the patient specimen by mixing **5 µl** of the patient specimen with **0.5 ml** of Serum Diluent.
5. Add **100 µl** of Calibrators, Positive and Negative controls and diluted patient specimens to the appropriate microwells indicated on the protocol record.

**Note:** Include one well with **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank. The absorbance of this well should not be greater than 0.3.

6. Incubate **30 minutes** (± 5 minutes) at room temperature on a level surface.

7. Wash step: Thoroughly aspirate the contents of each well. Add 200-300µL of the **reconstituted** wash buffer to all wells then aspirate. Repeat this sequence thrice more for a total of four washes. Invert the plate and tap it on absorbent material to remove any residual fluid after the last wash. Do not dry wells completely.
8. Add 100µL of the Conjugate to each well.
9. Incubate the wells for 30 minutes.
10. Wash step: Repeat step 7.
11. Add 100µL of Enzyme Substrate to each well.
12. Incubate for 30 minutes at room temperature.
13. Add 100µL of Stop Solution to each well. Maintain the same sequence and timing of Stop Solution addition as was used for the Enzyme Substrate. Read the absorbance (OD) of each well at 405nm within one hour of stopping the reaction.
14. Read the absorbance (OD) of each well at 405nm using a single or dual wavelength microplate reader against the reagent blank set at zero absorbance.

**Quality Control**

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <7 IU/ml. If the test is run in duplicate, take the mean of the two readings to determine the IU/ml. While performing Qualitative determinations, the optical density of the Calibrator D must be greater than that of the negative control and lesser than the absorbance of the positive control. For semi-quantitative determinations, the positive control must give values in the range stated on the vial.

**RESULTS**

**Calculations**

The concentrations of the patient samples can be determined by either of two methods:

**1. QUALITATIVE DETERMINATION**

**Abs. of Test Sample**

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{IU/ml of Calibrator D} = \text{IU/ml Test Sample}$$

**Abs. of Calibrator D**

**2. SEMI-QUANTITATIVE DETERMINATION**

Plot absorbance of Calibrator A through D against their respective concentration on a linear-linear graph paper. Plot the concentration in IU/ml on the X-axis against the absorbance on the Y-axis and draw the best fit curve. Determine the concentrations of the patient samples from the curve against its corresponding absorbance value.

**Calibrator**

The calibrators are included to provide semi-quantitation and must be used with each test. Patient specimens containing higher antibody levels may give absorbance values greater than that of the Calibrator A. For determining accurate semi-quantitative values such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining IU/ml, multiply the units obtained by the dilution factor.

**Interpretation**

The following information serves only as a guide in the interpretation of the laboratory results. Each laboratory must determine its own normal values.

RF-IgM	Interpretation/Interprétation/Interpretazione/
Value/Valueur/Valore/Valor/Wert	Interpretación/Deutung/Interpretação

<7 IU/ml	Neg (-)
7-9 IU/ml	Borderline/Indéterminé/Indeterminato/Indeterminado/Unbestimmt
>9 IU/ml	Pos (+)

The cutoff values for negative, indeterminate and positive were determined by testing 99 samples from healthy blood donors. When indeterminate values are obtained, the sample should be retested before reporting the results. If the values remain indeterminate, repeated testing of the patient ought to be considered.

#### **LIMITATIONS OF THE PROCEDURE**

The IgM-RF Test System should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only. The results obtained serve only as an aid in the diagnosis and should not be interpreted as diagnostic in themselves.

#### **EXPECTED VALUES**

RF is present in 59% of patients with rheumatoid arthritis (5). Serum IgM-RF is positive in 58% of Sjögren's Syndrome patients (8,9). As compiled from the cited literature, the incidence of RF in rheumatoid arthritis or other rheumatological disorders is summarized in table 1 at the end of this document:

#### **PERFORMANCE CHARACTERISTICS**

The IgM-RF Test System was compared with a commercially available test kit for the detection of RF in human serum.

A total of 87 sera were obtained from a clinical reference laboratory. They were identified as positive or negative for RF and were tested according to the procedures recommended by the manufacturer. The results are summarized in table 2 at the end of this document.

The ImmuLisa IgM-RF Test System showed a significant correlation and relative agreement with another predicate device, with a correlation coefficient of 0.8 and a slope of 0.95.

#### **Precision:**

The inter-assay and intra-assay Coefficient of Variation (CV) of the IgM-RF Test System ranges from 5% to 10% depending upon the levels of RF. Using three positive serum samples with varying RF levels in six replicates, intra- and inter- assay CV values are 0.8 to 4.1% and 3.0 to 8.2% respectively.

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## REFERENCES•REFERENCIAS•LITERATUR•RIFERIMENTI

1. Jonnson Thorbjorn, Valdimarsson H. Is measurement of rheumatoid factor isotypesclinically useful? Ann Rheum Dis 52: 161-164, 1993.
2. van Zeven D, Hazes JMW, Zwinderman AH, Cats A, Van der Voort EAM, Breedveld FC. Clinical significance of rheumatoid factors in earlyrheumatoid arthritis: results of a follow up study. Ann Rheum Dis 51:1029-1035, 1992.
3. Shmerling RH, Delbanco TL. The rheumatoid factor: an analysis of clinical utility. Am J Med 91:528-534, 1991.
4. Wolfe F, Cathey MA, Roberts FK. The latex test revisited-rheumatoid factor testing in 8,287 rheumatic disease patients. Arthr Rheum 34:951-960, 1991.
5. Bampton JLM, Cawston TE, Kyle MV, Hazelman BL. Measurement of rheumatoid factors by enzyme linked immunosorbant assay (ELISA) and comparison with other methods. Ann Rheum Dis 44: 13-19, 1985.
6. Markusse HM, Otten HG, Vroom TM, Smeets TJM, Fokkens N, Breedveld FC. Rheumatoid factor isotypes in serum and salivary fluid of patients with primary Sjögren's Syndrome. Ciin Immun Immunopath 66:26-32, 1993.
7. Atkinson JC, Travis WD, Slocum L, Ebbs WL, Fox P. Serum anti-SS-B/La and IgA rheumatoid factor are markers of salivary gland disease activity in primary Sjögren's syndrome. Arth Rheum 35:1368-1372, 1992.
8. Anderson LG, Talal N. The spectrum of benign to malignant lymphoproliferation in Sjögren's Syndrome. Clin Exp Immunol 10:199-221, 1972.
9. Saulsbury FT. Heavy and light chain composition of serum IgA and IgA rheumatoid factor in Henoch-Schönlein purpura Arthr Rheum 35:1377-1380, 1992.
10. Saulsbury FT. IgA rheumatoid factor in Henoch Schönlein purpura. J Pediatr 108:71-76, 1986.
11. Arnett FC. Revised criteria for the classification of rheumatoid arthritis. Bull Rheum Dis 38(5):1-6, 1989.
12. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, (HHS Pub No[CDC] 93-8395), 1993.

**Table 1. Incidence of RF in Patients with Rheumatic Diseases**

Rheumatoid Arthritis	70-90%
Sjögren's Syndrome	58%
Systemic lupus erythematosus	18%
Dermatomyositis	16%
Mixed connective tissue disease	13%
Cranial arteritis	10%
Polymyalgia rheumatica	10%
Polyarthritis	7%
Scleroderma	6%
Juvenile rheumatoid arthritis	6%

**Table 2. Comparison of Immulisa™ RF-IgM to another ELISA**

		IMMCO RF-IgM		
		POS (+)	NEG (-)	TOT (=)
<b>Other/Autre/Otro</b>	POS (+)	45	1	46
<b>Andere/Altro/Otra</b>	NEG (-)	2	39	41
<b>ELISA</b>	TOT (=)	47	40	87

relative specificity/spécificité relative/especificidad relativa/relative Spezifität/specificità relativa/especificidade relativa: 97%

relative sensitivity/sensibilitè relative/sensibilidad relativa/relative Sensitivität/sensibilità relativa/sensibilidade relativa: 98%

relative agreement/concordance/concordancia relativa/relative Übereinstimmung/concordanza relativa/acordo relativo: 95%

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**

**60 Pineview Drive**

**Buffalo, NY 14228-2120**

**Telephone: (716) 691-0091**

**Fax: (716) 691-0466**

**Toll Free USA/Canada: 1-800-537-TEST**

**E-Mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,

The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

[www.emergogroup.com](http://www.emergogroup.com)

REV.DEC2003

Document No. PI4139