# ImmuLisa™ Anti-human Tissue Tranglutaminase Antibody (hu tTG) IgA ELISA

For *in vitro* diagnostic use **IVD**

CLIA Complexity: High  
CDC Analyte Identification Code: 0546  
CDC Test System Identification Code: 28569

## PRODUCT INSERT

Catalog No. 1144

96 Determinations

## INTENDED USE

An enzyme linked immunosorbent assay (ELISA) for the detection and semi-quantitation of anti-human Tissue Transglutaminase IgG antibodies in human serum to aid in the diagnosis of IgA deficient patients with celiac disease.

## SUMMARY AND EXPLANATION

Celiac Disease (CD) is an autoimmune gastrointestinal disorder that may occur in genetically susceptible individuals triggered by the ingestion of gluten-containing grains such as wheat, barley and rye. Of the many autoimmune disorders, CD represents one of the few where the etiological agent is known and the disease subsides and goes in remission once the etiological agent is withdrawn from the diet.

CD is characterized by malabsorption resulting from inflammatory injury to the small intestinal mucosa and, when prolonged, can cause malnutrition. The classical symptoms of CD include diarrhea, weight loss and malnutrition.

Only a small percentage of patients with CD presents with classical symptoms. Consequently, the clinical spectrum of CD has grown much broader than in the past to include patients that do not present with classical symptoms. It is not uncommon for the initial symptoms to be non-gastrointestinal or for gastrointestinal symptoms, if present, to be mild or intermittent.

Some common non-gastrointestinal manifestations include short stature, iron and folate deficiency, anemia, bone loss, aphthous stomatitis, arthralgia, dental enamel defects, etc. The need to examine a wider range of clinical presentation has led to greater numbers of individuals diagnosed with CD later in life than ever before. Adults may present with iron deficiency, macrocytic anemia and hypocalcaemia.

Studies have found the prevalence of CD to be highly variable from population to population.<sup>1</sup> The true prevalence has been difficult to ascertain. Disparate criteria used in diagnosis of CD are often the cause. If clinical criteria alone are used in determining prevalence, the incidence of CD is much lower as compared with incidence established by serological methods.<sup>1,2</sup> Using serological methods, the incidence of CD in the general population is approximately one in 200.

Diagnosis of CD based on clinical criteria thus can be misleading and may lead to serious delays in proper diagnosis. Frequently, delays in diagnosis extend 10-13 years from the first clinical presentation of symptoms.

Failure to diagnose CD early on may predispose an individual to long-term complications such as splenic atrophy and intestinal lymphoma. The incidence of lymphoma arising in the context of CD is difficult to ascertain. One study has shown incidence of lymphoma involving the gastrointestinal (GI) tract in patients with CD to range from 3.6 percent to 40 percent.<sup>3</sup> In another recent study,<sup>4</sup> CD is associated with significantly elevated risk for intestinal

lymphoma, especially for non-Hodgkin's. A gluten-free diet (GFD) normalizes the mucosa and helps reduce the malignant potential.

Histological examination of the small intestinal biopsy remains the gold standard for diagnosing CD, but it has its own limitations. These include some patients with latent or even active CD that may have normal histopathology.<sup>5</sup>

The revised European Society of Paediatric Gastroenterology and Nutrition (ESPGHAN) criteria include only a single biopsy with clear-cut remission of clinical symptoms on GFD<sup>6</sup>. Positive serology at the time of diagnosis with disappearance on GFD contributes to the diagnosis. The various serological tests employed in the work-up of patients suspected to have CD include anti-gliadin antibody (AGA), anti-endomysial antibody (EMA), anti-reticulin antibody (ARA) and anti-tissue transglutaminase (tTG) antibody tests. Antibodies to gliadin and tTG are detected by ELISA, whereas endomysium and reticulin antibodies are detected by indirect immunofluorescence. EMA are very specific indicators of CD. However, the EMA test is an immunohistochemical method that requires experience in reading immunofluorescence reactions. Rosario et al<sup>7</sup> concludes that "EMA-IgA is 100 percent sensitive and specific in active, untreated IgA-sufficient CD patients when performed by an established laboratory."

Since identification of tTG as the endomysial antigen, ELISA methods have been described for detecting antibodies in the sera of patients with CD. The advantage of the anti-tTG antibody assay is that it is automatable and less subjective than EMA. For this reason, many laboratories have opted to use the tTG antibody method as the screening method. In these laboratories, it may be the only assay used for detection of CD cases. In various studies on the efficacy of the tTG antibody method for screening for CD, the specificity and sensitivity of this method has been found to range from 90 percent to 95 percent.<sup>8</sup> Human tTG has been described to improve the sensitivity of the tTG antibody assay for CD.

## PRINCIPLES OF PROCEDURES

Recombinant human Tissue Transglutaminase antigen is bound to the wells of a polystyrene microwell plate followed by blocking the unreacted sites to reduce non-specific binding. Controls, calibrators and diluted patient sera are added to separate wells, allowing any anti-tTG antibodies present to bind to the immobilized antigen. Unbound sample is washed away and an enzyme labeled anti-human IgA conjugate is added to each well. These enzyme conjugated antibodies bind specifically to the human immunoglobulin of the IgA class. After washing away any unbound conjugate, specific enzyme substrate (pNPP) is then added to the wells. After stopping the enzymatic reaction, the intensity of color change, which is proportional to the concentration of antibody, is read by a spectrophotometer at 405 nm. Results are expressed in ELISA units per milliliter (EU/ml).

## REAGENTS

### Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Do not use if reagent is not clear or if a precipitate is present. All reagents must be brought to room temperature (20-25°C) prior to use. When stored at 2-8°C, the reconstituted wash buffer is stable until the kit expiration date. Reconstitute the wash buffer to 1 liter with distilled or deionized water. Coated microwell strips are for one time use only.

### Precautions

For *in vitro* Diagnostic Use. All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>9</sup>. **WARNING** - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

**Instructions should be followed exactly as they appear in this kit insert to ensure valid results.** Do not interchange kit components with those from other sources other than the same catalog number from IMMCO DIAGNOSTICS.

Use good laboratory techniques to minimize microbial and chemical contamination.

Do not use after expiration date.

**Materials Provided**

ImmuLisa™ IgA-tTG ELISA REF 1144

Kits contain sufficient reagents to perform 96 determinations each.

12 x 8	<b>MICROPLATE</b>   tTG	<b>Microplate</b> with individual breakaway microwells coated with hu tTG antigen.
1 x 1.5 ml	<b>CALIBRATOR A</b>   tTG * †	Ready to use <b>Calibrator A</b> ( <i>green cap</i> ). Human serum containing antibodies to tTG.
1 x 1.5 ml	<b>CALIBRATOR B</b>   tTG * †	Ready to use <b>Calibrator B</b> ( <i>violet cap</i> ). Human serum containing antibodies to tTG.
1 x 1.5 ml	<b>CALIBRATOR C</b>   tTG * †	Ready to use <b>Calibrator C</b> ( <i>blue cap</i> ). Human serum containing antibodies to tTG.
1 x 1.5 ml	<b>CALIBRATOR D</b>   tTG * †	Ready to use <b>Calibrator D</b> ( <i>yellow cap</i> ). Human serum containing antibodies to tTG.
1 x 1.5 ml	<b>CONTROL +</b>   tTG * †	Ready to use <b>Positive Control</b> ( <i>red cap</i> ). Contains human serum positive for tTG.
1 x 1.5 ml	<b>CONTROL -</b> *	Ready to use <b>Negative Control</b> ( <i>white cap</i> ). Contains human serum.
1 x 12 ml	<b>IgA-CONJ</b>   <b>ALKPHOS</b> * †	Ready to use <b>anti-human Alk. Phos. Conjugate</b> . Color coded pink.
1 x 60 ml	<b>DIL</b> *	Ready to use <b>Serum Diluent</b> . Color coded blue.
1 x 12 ml	<b>SUBSTRATE</b> *	Ready to use <b>Enzyme Substrate</b> . Contains pNPP. <b>Protect from light.</b>
1 x 12 ml	<b>STOP</b>	Ready to use <b>Stop Solution</b> .
2	<b>BUF</b>   <b>WASH</b>	Powder <b>Wash Buffer</b> . Reconstitute to one liter each.

\* Contains < 0.1% NaN<sub>3</sub>

**Materials Required But Not Provided**

- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Deionized or distilled water
- Microplate reader capable of reading absorbance values at 405 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600-650 nm.
- Squeeze bottle to hold diluted wash buffer
- Timer
- Absorbent paper
- Automatic microplate washer capable of dispensing 200 µl

**SPECIMEN COLLECTION AND HANDLING**

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°- 8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

## PROCEDURE

### Procedural Notes

- Before starting the assay read these instructions carefully.
- Bring all reagents and samples to room temperature (20-26°C) for 30 minutes. Return materials to refrigerator immediately after use.
- Prepare all dilutions of the patient specimens before starting the test.
- **Immediately return unused strips to the pouch containing desiccants and seal securely to minimize exposure to water vapor.**
- Wash step: Good technique is critical. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 wells simultaneously. This speeds the process and provides for a more uniform incubation time.
- Careful timing is important. Incubation periods begin after dispensing reagents.

### Assay procedure

1. **ALL REAGENTS MUST BE BROUGHT TO ROOM TEMPERATURE (20-26°C) PRIOR TO BEGINNING THE ASSAY.**
2. Label protocol record to indicate specimen placement in the microplate. It is good laboratory practice to test specimens in duplicate.
3. **Qualitative determination:** use only Calibrator D.  
**Semi-quantitative determination:** use Calibrators A - D as shown in the example below.

#### QUALITATIVE DETERMINATION

A	BLANK	S5		
B	NEG	S6		
C	POS	S7		
D	CAL D	S8		
E	S1	S9		
F	S2	S10		
G	S3	S11		
H	S4	S12		
	1	2	3	4

#### SEMI-QUANTITATIVE DETERMINATION

A	BLANK	S2		
B	NEG	S3		
C	POS	S4		
D	CAL A	S5		
E	CAL B	S6		
F	CAL C	S7		
G	CAL D	S8		
H	S1	S9		
	1	2	3	4

4. Prepare a **1:51** dilution of the patient specimen by mixing **10 µl** of the patient specimen with **0.5 ml** of Serum Diluent.
5. Add **100 µl** of Calibrators, Positive and Negative controls and diluted patient specimens to the appropriate microwells indicated on the protocol record.  
**Note:** Include one well with **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank. The absorbance of this well should not be greater than 0.3.
6. Incubate **30 minutes** ( $\pm$  5 minutes) at room temperature on a level surface.
7. Wash step: Thoroughly aspirate the contents of each well. Add 200-300µL of the **reconstituted** wash buffer to all wells then aspirate. Repeat this sequence thrice more for a total of four washes. Invert the plate and tap it on absorbent material to remove any residual fluid after the last wash. Do not dry wells completely.
8. Add 100µL of the Conjugate to each well.
9. Incubate the wells for 30 minutes.
10. Wash step: Repeat step 7.
11. Add 100µL of Enzyme Substrate to each well.

12. Incubate for 30 minutes at room temperature.
13. Add 100µL of Stop Solution to each well. Maintain the same sequence and timing of Stop Solution addition as was used for the Enzyme Substrate. Read the absorbance (OD) of each well at 405nm within one hour of stopping the reaction.
14. Read the absorbance (OD) of each well at 405nm using a single or dual wavelength microplate reader against the reagent blank set at zero absorbance.

**Quality Control**

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. Calibrator A should have an absorbance reading  $\geq 1.0$ , otherwise the test must be repeated. The negative control must be <20 EU/ml. If the test is run in duplicate, take the mean of the two readings to determine the concentration of anti-hu tTG antibodies. While performing qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the positive control. For semi-quantitative determinations, the positive control must give values in the range stated on the vial. We recommend borderline samples be tested with a fresh sample taken at a later date to ensure accuracy.

**RESULTS**

**Calculations**

The concentrations of the patient samples can be determined by either of two methods:

**1. QUALITATIVE DETERMINATION**

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

**2. SEMI-QUANTITATIVE DETERMINATION**

Plot absorbance of Calibrator A through D against their respective concentration on a linear-linear graph paper. Plot the concentration in EU/ml on the X-axis against the absorbance on the Y-axis and draw the best fit curve. Determine the concentrations of the patient samples from the curve against its corresponding absorbance value. See Figure 1 at the end of this document for a sample standard curve.

**Calibrator**

The calibrators are included to provide semi-quantitation and must be used with each test. Patient specimens containing higher antibody levels may give absorbance values greater than that of the Calibrator A. For determining accurate semi-quantitative values such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml, multiply the units obtained by the dilution factor.

**Interpretation**

The following information serves only as a guide in the interpretation of the laboratory results. Each laboratory must determine its own normal values.

Anti-tTG Value/Valeur/Valore/Valor/Wert	Interpretation/Interprétation/Interpretazione/ Interpretación/Deutung/Interpretação
<20 EU/ml	Neg (-)
20 – 25 EU/ml	Borderline/Indéterminé/Indeterminato/ Indeterminado/Unbestimmt
>25 EU/ml	Pos (+)

**LIMITATIONS OF THE PROCEDURE**

The Immulisa™ anti-tTG IgA assay should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. The method should be used for testing human serum samples only. Results obtained serve only as an aid in the diagnosis and should not be interpreted as

diagnostic in themselves. In IgA deficient CD patients anti-tTG may be negative.

### **EXPECTED VALUES**

The expected values in a normal population are negative (<20 EU/ml for adults and children). However, it has been determined that some apparently healthy, asymptomatic individuals may test positive for IgA anti-tTG antibodies.

The incidence of anti-tTG antibodies has been studied by various investigators and the findings are summarized in Table 1 at the end of this document.

The incidence and the levels of anti-tTG antibodies is dependent upon the diet status. The levels of these antibodies decrease and eventually will become negative in patients with CD who are on a gluten-free diet. Similarly, the levels of these antibodies will increase or the anti-tTG antibodies become positive when patients with CD who were on a gluten-free diet ingests a gluten-containing diet <sup>10, 11</sup>.

As there are several antibody markers associated with CD, the algorithm in Figure 2 at the end of this document may serve as a guide in the interpretation of various serological markers for CD.

Some investigators have taken the approach presented in Figure 3 to establish a diagnosis of CD.

### **PERFORMANCE CHARACTERISTICS**

The utility of the ImmuLISA™ anti-tTG IgA Antibody ELISA was determined by comparing the results with:

- a) another commercially available anti-tTG IgA ELISA method and
- b) ImmuGlo™ anti-endomysial antibody immunofluorescence method.

Normal Range: The normal range was established by testing 64 serum samples from apparently healthy donors obtained from the Red Cross. The cutoff equals the average value + 3 STD of these samples.

### **Comparative Specificity and Sensitivity**

- A. Endomysial Antibodies vs. ImmuLISA™ tTG IgA Antibody ELISA: A total of 274 samples were tested for EMA using ImmuGlo™ EMA kit and the results were compared with ImmuLISA anti-tTG IgA ELISA method. The results are summarized in Table 2:
- B. ImmuLISA™ anti-tTG ELISA vs. Another Commercial anti-tTG IgA ELISA Method: A total of 74 samples were tested on ImmuLISA™ anti-tTG IgA antibody kit and another commercially available anti-tTG IgA kit. The results appear in Table 3:
- C. Cross Reactivity: A total of 65 disease controls from autoimmune vesiculo-bullous cases such as pemphigus were tested. None were found positive.

### **Precision:**

Samples with known concentrations of anti-tTG IgA were assayed in 10 replicates over a period of two weeks. Intra-and inter-assay coefficient of variation (CV) were calculated. See Table 4.

### **Recovery:**

Samples with known anti-tTG IgA concentrations were mixed with appropriate dilutions of another positive sample with known amounts of anti-tTG IgA. Anti-tTG IgA antibody levels of the mixed samples were determined and from the values obtained the percent recovery calculated. The results appear in Table 5.

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
## REFERENCES•REFERENCIAS•LITERATUR•RIFERIMENTI

1. Cooper BT, Holmes GK, Ferguson R, et al. Celiac disease and malignancy. *Medicine*; 1980, 59:249-261.
2. Sher KS, Jayanthi V, Probert CS, et al. Infertility, obstetric and gynecological problems in coeliac sprue. *Digestive Dis*; 1994, 12:186-190.
3. McFarlane XA, Bhalla AK, Reeves DE, et al. Osteoporosis in treated adult coeliac disease. *Gut*; 1995, 36:710-714.
4. Gobbi G, Bouquet F, Greco L, et al. Coeliac disease, epilepsy and cerebral calcifications. The Italian working group on coeliac disease and epilepsy. *Lancet*; 1992, 340:439-443.
5. Dietrich W, Ehnis T, Bauer M, et al. Identification of tissue transglutaminase as the autoantigen of celiac disease. *Nature Medicine*; 1997, 3:797-801.
6. Phillips MA, Stewart BE, Quin Q, et al. Primary structure of keratinocyte transglutaminase. *Proc Natl Acad Sci*; 1990, 87:9333-9337.
7. Aschlimann D, Mosher D, Paulsson M, et al. Tissue transglutaminase and factor XIII in cartilage and bone remodelling. *Sem Thromb and Hemo*; 1996, 22:437-443.
8. Jones RA, Nicholas B, Mian S, et al. Reduced expression of tissue transglutaminase in a human endothelial cell line leads to changes in cell spreading, cell adhesion and reduced polymerization of fibronectin. *J Cell Sci*; 1997, 110:2461-2472.
9. Sulkanen S, Haltunen T, Laurila K, et al. Tissue transglutaminase autoantibody enzyme-linked immunosorbent assay in detecting celiac disease. *Gastroenterology*; 1998, 115:1322-1328.
10. Dietrich W, Laag E, Schopper H, et al. Autoantibodies to tissue transglutaminase as predictors of celiac disease; *Gastroenterology* 1998, 115:1317-1321.
11. Molberg O, McAdam S, Korner R. Tissue transglutaminase selectively modifies gliadin peptides that are recognized by gut-derived T-cells in celiac disease. *Nature Medicine*; 1998, 4:713-717.
12. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health; 1993 (HHS Pub No [CDC] 93-8395).
13. Dietrich W, Laag E, Bruckner-Tuderman et al. Antibodies to tissue transglutaminase as serologic markers in patients with dermatitis herpetiformis. *J Invest Dermatol*; 1999, 113:133-136.

## Figure 1. Sample Standard Curve

Standard Curve/Courbe d'étalonnage/Curva standard/  
Curva estándar/Standardkurve/Curva padrão  
**Anti-hu tTG Immulisa™**

**Table 1. Incidence of EMA and tTG Antibodies**


Study Group	No.	EMA	tTG
Dermatitis herpetiformis	44		
on gluten diet	33	33 (100%)	33 (100%)
on gluten-restricted diet	11	0 (0%)	3 (27%)
Celiac Disease	58		
on gluten diet	34	34 (100%)	31 (91%)
on gluten-restricted diet	24	0 (0%)	1 (4%)
Controls	189	0 (0%)	0 (0%)
Pemphigus	80	0 (0%)	1 (1%)
Linear IgA dermatosis	3	0 (0%)	0 (0%)
Normal	106	0 (0%)	1 (1%)

*Kumar V et al; Clin Immunol 2001; 98: 378-382*

**Figure 2: CD Diagnosis Algorithm**

EMA	ARA	AGA		anti-tTG	Interpretation
		IgG	IgA		
+	+	+	+	+	Positive
+	+	-	+	+	Positive
+	+	-	-	+	Positive
+	-	+	+	+	Positive
+	-	-	-	+	CD – probable
-	+	+	+	+	CD – probable
-	-	+	+	+	CD – possible
-	-	+	-	-	CD – unlikely
-	-	-	+	-	CD – unlikely
-	-	-	-	-	Negative

**Figure 3: Laboratory Testing Algorithm**

**Table 2. ImmuLisa™ anti-hu tTG IgA vs ImmuGlo™ EMA**

		anti-hu tTG IgA ELISA		
		Pos (+)	Neg (-)	Total (=)
ImmuGlo™ EMA	Pos (+)	109	3	112
	Neg (-)	5	157	162
	Total (=)	114	160	274

relative specificity/spécificité relative/especificidad relativa/relative Spezifität/specificità relativa/especificidade relativa: 97%  
 relative sensitivity/sensibilit  relative/sensibilidad relativa/relative Sensitivit t/sensibilit  relativa/sensibilidade relativa: 98%  
 relative agreement/concordance/concordancia relativa/relative  bereinstimmung/concordanza relativa/acordo relativo: 94%

**Table 3. ImmuLisa™ anti-hu tTG IgA vs other anti-hu tTG IgA ELISA**

		ImmuLisa™ anti-hu tTG IgA ELISA		
		Pos (+)	Neg (-)	Total (=)
<b>Other/Autre/Otro</b>	Pos (+)	57	1	58
<b>Anderes/Altri/Outro</b>	Neg (-)	1	15	16
<b>ELISA</b>	Total (=)	58	16	74

relative specificity/spécificité relative/especificidad relativa/relative Spezifität/specificità relativa/especificidade relativa: 97%  
 relative sensitivity/sensibilità relative/sensibilidad relativa/relative Sensitivität/sensibilità relativa/sensibilidade relativa: 96%  
 relative agreement/concordance/concordancia relativa/relative Übereinstimmung/concordanza relativa/acordo relativo: 98%

**Table 4. Precision**

Anti-hu tTG IgA	Intra-plate	Inter-plate
	plaque/piastra/placa/Platte CV	plaque/piastra/placa/Platte CV
1. 199 EU/ml	6.2%	2.9%
2. 85 EU/ml	8.4%	6.5%
3. 41 EU/ml	15.6%	9.6%

**Table 5. Recovery**

	EU/ml added/supplémentaire/ aggiunto/agregado/ hinzugefügt/adicionado	EU/ml obtained/obtenuto/ verificato/obtenido/ erreicht/obtido	% Recovery/Rétablissement/ Recupero/Recuperación/ Wiederaufnahme/Recuperação
1.	107.5	112	104.2%
2.	86	94.6	110%
3.	64	70.9	109.9%

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**  
**60 Pineview Drive**  
**Buffalo, NY 14228-2120**  
**Telephone: (716) 691-0091**  
**Fax: (716) 691-0466**  
**Toll Free USA/Canada: 1-800-537-TEST**  
**E-Mail: info@immcodiagnostics.com**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
 Autorizzato/Representante Autorizado/Représentant Autorisé  
 EMERGO Group, Inc.  
 Molenstraat 15, 2513 BH, The Hague,  
 The Netherlands  
 Tel (+31) 345 8570, Fax (+31) 346 7299  
 www.emergogroup.com