ImmuGlo™ Anti-Endomysial Antibody (EMA) Test System

For *in vitro* diagnostic use

IVD

CLIA Complexity: High
CDC Analyte Identification Code: 0497
CDC Test System Identification Code: 28281

PRODUCT INSERT

Code 1114
Code 1114-96

48 Determinations
96 Determinations

INTENDED USE

An indirect immunofluorescence antibody test for the qualitative and semi-quantitative detection of endomysial antibodies (EMA) in human serum as an aid in the diagnosis of celiac disease and dermatitis herpetiformis.

SUMMARY AND EXPLANATION

Endomysial antibodies (EMA), as reported in the literature, are detected primarily on the smooth muscle of monkey esophagus by indirect immunofluorescence. The detection of EMA aids in the diagnosis of *gluten sensitive enteropathy*, i.e. *celiac disease* (CD) and *dermatitis herpetiformis* (DH). Patients with CD and DH are reported to have antibodies to endomysium, reticulin and gliadin¹⁻¹². These serological markers have recently been incorporated into the revised criteria for the diagnosis of CD by the European Society of Pediatric Gastroenterology and Nutrition¹³. Of the various antibody markers of CD and DH, EMA of the IgA class seem to be the most sensitive and specific marker. EMA of the IgG class also occur when IgA class EMA are in high titer or in individuals who are IgA deficient. A rapid decrease in EMA levels results with adherence to a gluten free diet. A gluten challenge or a failure to maintain a gluten free diet leads to the appearance or an increase in endomysial antibody titers. Patients on a gluten free diet >9 months have reduced or negative EMA titers if they adhere to their diet restrictions^{1,6-8,10}.

PRINCIPLES OF PROCEDURES

In the indirect immunofluorescence method, patient serum is incubated on tissue sections to allow binding of antibodies to the substrate. Any antibodies not bound are removed by rinsing. Bound antibodies of the IgA and IgG class are detected by incubation of the substrate with fluorescein-labeled, anti-human immunoglobulin conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of EMA is demonstrated by an apple green fluorescence of the endomysial lining of smooth muscle bundles. The titer (the reciprocal of the highest dilution giving a positive reaction) of the antibody is then determined by testing serial dilutions¹⁴

PRODUCT INFORMATION

Storage and Preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

Materials Provided

ImmuGlo™ EMA IFA REF 1114, 1114-96

Kits contain sufficient reagents to perform 48 determinations each.

8 x	SORB SLD	6 well Primate Smooth Muscle Substrate Slides
16 x	SORB SLD	6 well Primate Smooth Muscle Substrate Slides (1114-96)
1 x 0.5 ml	CONTROL + EMA *	EMA Positive Control. Contains human serum.
1 x 0.5 ml	CONTROL - *	Negative Control. Contains human serum.
1 x 5 ml	IgA/IgG-CONJ FITC *	Anti-human polyvalent FITC Conjugate. Protect from light. 1 vial (1114), 2 vials (1114-96)
1 x 60 ml	BUF *	Sample Diluent.
2 vials	BUF WASH	Phosphate Buffered Saline (PBS). Dissolve each vial to 1 liter.
1 x 5.0 ml	MOUNTING MEDIUM *	Mounting Medium. Do not freeze. 1 vial (1114), 2 vials (1114-96)
1 x 1.0 ml	EVANS	Evan's Blue Counterstain.
1 x 12	COVER SLD	Coverslips 1 box (1114), 2 boxes (1114-96)

* Contains < 0.1% NaN₃

Materials Required But Not Provided

Fluorescence microscope
Micropipette or Pasteur pipette
Serological pipettes
Staining dish (e.g. Coplin jar)
Small test tubes (e.g. 13 x 75 mm) and test tube rack
Distilled or deionized water
1 liter container
Wash bottle
Paper towels
Incubation chamber

WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials¹⁴.

WARNING - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from IMMCO. Do not use beyond expiration date.

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

PROCEDURE

Test Method

A. Screening

1. Dilute each patient serum 1:2.5 with the Sample Diluent provided (0.2 ml serum + 0.3 ml Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 50 µl of the Negative Control to well #1. Similarly apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 50 µl) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 50 µl) to each well.
9. Repeat steps **7 and 8** for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If desired, 2-3 drops of Evans blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.
12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

B. Endpoint Determination (titration)

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Each test run should include the Positive and Negative Controls. Make serial two-fold dilutions starting at 1:2.5. The reciprocal of the highest dilution producing a positive reaction is the titer.

Preparation of Serial Dilutions

Number four tubes 1 through 4. Add 0.3 ml of Buffered Diluent to tube 1 and 0.2 ml to tubes 2 through 4. Pipette 0.2 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

Tubes	1	2	3	4
Serum	0.2 ml			
	+			
Buffered Diluent	0.3 ml	0.2 ml	0.2 ml	0.2 ml
Transfer		0.2 ml	0.2 ml	0.2 ml
Final dilution	1:2.5	1:5	1:10	1:20 etc.

QUALITY CONTROL

Both a Positive and Negative Control should be included with each test run. The Negative Control should show no specific fluorescence of the endomysium lining of the smooth muscle bundles, whereas the Positive Control should have 2+ or greater staining intensity of the tubules of these structures.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control
- Problems with the optical system of the fluorescence microscope. These may include : improper alignment, bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

RESULTS

The results of the tests for endomysial antibodies should be reported as negative (<2.5), positive greater or equal to 20, or preferably, positive with titer.

Read for specific staining of the endomysium lining of the smooth muscle bundles. **See Photo 1.** Endomysial antibodies react as a network of thin, irregular lines around the sarcolemma of the individual smooth muscle fibrils. This is in a sharp contrast to anti-smooth muscle antibodies which react with the sarcoplasm. **See Photo 2.**

Other detectable antibodies besides anti-smooth muscle antibodies (ASMA) include anti-nuclear antibodies (ANA). The presence of ASMA is known to cause false negative results for endomysial antibodies. If ASMA are detected, then the sample should be tested at higher dilutions¹. ANA reactions on smooth muscle tissue, when they occur, are usually weak and sparsely distributed and, therefore, unlikely to cause false negative results.

Consult Photo 1 and Photo 2 at the end of this document for example reactions.

LIMITATIONS OF THE PROCEDURE

The ImmuGlo™ EMA IFA Test System: CLIA Complexity: High; CDC Analyte Identification Code: 0497; CDC Test System Identification Code: 2828.

In some cases, sera positive for EMA may either be very weak or negative at the initial screening dilution (prozone phenomenon). In such doubtful cases the sera should be screened at higher dilutions and, if positive, antibody titer determined.

The presence of two or more antibodies in a serum which are reactive with the same tissue may cause an interference in their detection by immunofluorescence. This interference may cause either a failure to detect EMA or a suppression of its titer if the interfering antibody has a higher titer than EMA. The most common cause of the interference phenomenon in EMA tests is the coexistence of smooth muscle antibodies. It is recommended that patients sera which also contain ASMA be tested further at higher dilutions. IgA class ASMA are not a common occurrence. IgG class ASMA do not block IgA-EMA as the former react with the sarcoplasm of smooth muscle bundles and the latter react with the endomysium of the sarcolemma around the smooth muscle bundles. Anti-reticulin antibodies do not interfere with the reaction of EMA because they do not react with primate smooth muscle tissue. The coexistence of IgG class EMA may interfere with the detection of IgA class EMA. However, this rarely occurs as:

1. IgG class - EMA are present in only 25% of celiac disease patients,
2. IgG class - EMA titers are usually much lower than IgA-EMA titers and
3. IgA antibodies are usually of higher avidity than IgG antibodies.

In some patients with celiac disease and IgA deficiency, the IgA-class endomysial antibodies are absent. However, such patients are usually positive for IgG class EMA.

Patients with celiac disease on a gluten free diet for >9 months invariably are negative for EMA.

When making a diagnosis, results of all laboratory testing must always be evaluated along with the total clinical history of the patient.

EXPECTED VALUES

As seen in Table 1, EMA, as detected on primate smooth muscle are highly specific markers for celiac disease and dermatitis herpetiformis. The presence of EMA seems to be related to the intestinal pathology both in celiac disease and dermatitis herpetiformis rather than to the skin lesions in the latter, as depicted in Figure 1.

PERFORMANCE CHARACTERISTICS

The ImmuGlo™ Anti-Endomysial Antibody (EMA) Test kit, using primate smooth muscle substrate and a polyvalent conjugate, was compared with another commercially available kit also using a polyvalent conjugate and monkey esophagus as a substrate. The comparison included a total of 68 sera: 20 from patients with clinically suspected celiac disease and 48 from normal controls. Sera were tested according to the procedure recommended by the manufacturer. A screening dilution of 2.5 was used and all sera positive for EMA were titrated to endpoint. The results were as follows:

		IMMCO EMA		
		POS (+)	NEG (-)	TOT (=)
Other/Autre/Otro	POS (+)	18	0	18
Andere/Altro/Otra	NEG (-)	2	48	50
EMA IFA	TOT (=)	20	48	68

relative specificity/spécificité relative/especificidad relativa/relative Spezifität/specificità relativa/especificidade relativa: 97%
relative sensitivity/sensibilit  relative/sensibilidad relativa/relative Sensitivit t/sensibilit  relativa/sensibilidade relativa: 100%
relative agreement/concordance/concordancia relativa/relative  bereinstimmung/concordanza relativa/acordo relativo: 96%

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8x	SORB SLD 6
16x	SORB SLD 6
1 x 0,5 ml	CONTROL + EMA *
1 x 0,5 ml	CONTROL - *
1 x 5 ml	IgA/IgG-CONJ FITC *
1 x 60 ml	BUF *
2 vials	BUF WASH
1 x 5,0 ml	MOUNTING MEDIUM *
1 x 1,0 ml	EVANS
1 x 12	COVER SLD

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REFERENCES • REFERENCIAS • LITERATUR • RIFERIMENTI

1. Chorzelski TP, Beutner EH, Kumar V and Zalewski TK (Eds). Serologic Diagnosis of Celiac Disease. CRC Press Boca Raton; 1990.
2. Levenson SD, Austin RK, Dietler MD, Kasarda DD and Kagnoff MF. Specificity of antigliadin antibody in celiac disease. *Gastroenterology*; 1985, 89:1-5.
3. Tucker NT, Barghuthy FS, Prihoda TJ, Kumar V, Kerner A and Lebenthal E. Antigliadin antibodies detected by enzyme linked immunosorbent assay as a marker of childhood celiac disease. *J Pediatric*; 1988, 113:286-289.
4. Hallström O. Comparison of IgA-class reticulin and endomysium antibodies in coeliac disease and dermatitis herpetiformis. *Gut*; 1989, 30:1225-1232.
5. Khoshoo V, Bhan MK, Unsworth DJ, Kumar R and Walker-Smith A. Anti-reticulin antibodies: useful adjunct to histopathology in diagnosing celiac disease, especially in developing country. *J Pediatric Gastroenterol Nutr*; 1988, 7: 864-866.
6. Kapuscinska A, Zalewski T, Chorzelski TP, Sulez J, Beutner EH, Kumar R and Rossi T. Disease specificity and dynamics of changes in IgA class anti-endomysial antibodies in celiac disease. *J Pediatric Gastroenterol Nutr*; 1987, 6: 529-534.
7. Rossi TM, Kumar V, Lerner A et al. Relationship of endomysial antibodies to jejunal mucosal pathology: specificity towards both symptomatic and asymptomatic celiacs. *J Pediatric Gastroenterol Nutr*; 1988, 7: 858-863.
8. Kumar V, Lerner A, Valeski JE et al. Endomysial antibodies in the diagnosis of celiac disease and the effects of gluten on antibody titers. *Immun Invest*; 1989 18: 533-544.
9. Calabuig M, Torregosa R, Polo P, Tuset L et al. Serological markers and celiac disease: a new diagnostic approach. *J Paediatric Gastroenterol Nutr*; 1990, 10:435-442.
10. Karpati S, Bürgin-Wolff A, Krieg T, Maurer M et al. Binding to human jejunum of serum IgA antibody from children with coeliac disease. *Lancet*; 1990, 336: 1335-1338.
11. Kumar V, Hemedinger E, Chorzelski TP et al. Reticulin and endomysial antibodies in bullous disease - comparison of specificity and sensitivity. *Arch Dermatol*; 1987, 123: 1179-1182.
12. Peters MS and McEvoy MT. IgA antiendomysial antibodies in dermatitis herpetiformis. *J Am Acad Dermatol*; 1989, 21: 1225-1231.
13. Walker-Smith JA, Quandalini S, Schmitz J et al. Revised criteria for diagnosis of celiac disease. Report of working group of European Society of Paediatric Gastroenterology and Nutrition. *Arch Dis Child*; 1990, 65:909-911.
14. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons New York; 1987, 3rd Ed, 3-40.
15. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1993 (HHS Pub. No. {CDC} 93-8395).

Photo 1. EMA staining reaction on primate smooth muscle, 200X. Note staining of lining of the smooth muscle bundles.


Photo 2. ASMA staining reaction on primate smooth muscle, 200X. Note staining of the smooth muscle sarcoplasm.



Figure 1: Correlation of IgA-EMA titers in Villous Atrophy

From Chorzelski TP et al¹ and Kumar V et al⁸

Table 1. Incidence of IgA Class EMA

Clinical Condition	No. Tested	% Positive
Confirmed Celiacs		
On gluten	185	99
On gluten free diet	190	9
Suspected Celiacs		
On gluten	82	83
On gluten free diet	30	16
Dermatitis Herpetiformis (DH)	253	80
DH with Subtotal Villous Atrophy	42	100
DH on gluten free diet	36	3
Disease Controls (GI)		
Infectious Diarrhea	210	0
Recurrent Diarrhea	124	0
Toddlers Diarrhea	170	0
Milk Protein Intol.	69	0
Ulcerative Colitis	69	0
Crohn's Disease	65	0
Liver Diseases	21	0
Disease Controls (Skin)		
Linear IgA Bullous Dermatitis	4	0
Other Skin Diseases	180	0

Compiled from the literature as per Chorzelski TP et al¹.

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante
 Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.
 Molenstraat 15, 2513 BH, The Hague,
 The Netherlands
 Tel (+31) 345 8570, Fax (+31) 346 7299
www.emergogroup.com

For technical assistance please contact:

IMMCO Diagnostics, Inc.
60 Pineview Drive
Buffalo, NY 14228-2120
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immcodiagnostics.com

or your local product distributor