# ImmunoGlo™ Anti-Thyroid Antibody (ATA) Test System

**IVD**

**REF** Code: 1143 48 determinations

## PRODUCT INSERT

### INTENDED USE

An indirect immunofluorescence antibody test for the detection and semi-quantitation of thyroid autoantibodies in human serum.

### SUMMARY AND EXPLANATION

Autoimmune diseases such as Hashimoto's thyroiditis, primary autoimmune hypothyroidism and Grave's disease, are characterized by the presence of autoantibodies to thyroglobulin and microsomal antigens. Thyroglobulin is a large (660 kD) glycoprotein and functions as a prohormone. The microsomal antigen, identified as thyroid peroxidase (TPO) 105 kD antigen, is involved in the iodination and coupling of specific homogenous tyrosines in the production of thyroxine and triiodothyronine. Measurement of autoantibodies to thyroglobulin and microsomal (TPO) antigens are important in the diagnosis of thyroid diseases<sup>1-4</sup>. The thyroid antibodies can be measured by various methods such as indirect immunofluorescence, passive hemagglutination, and enzyme immunoassays.

### PRINCIPLES OF PROCEDURE

In the indirect immunofluorescence method used, patient serum is incubated on optimized preparations of tissue sections to allow binding of antibodies to the substrate. Any antibodies not bound are removed by rinsing. Bound antibodies of the IgG class are detected by incubation of the substrate with fluorescein-labeled anti-human IgG conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters<sup>5</sup>.

The presence of autoantibody reactions are demonstrated by an apple green fluorescence, either of the colloid or the follicular epithelial cells.

### PRODUCT INFORMATION

#### Storage and preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

#### Materials provided

**REF** Code: 1143 48 determinations

8 x	<b>SORB SLD 6</b>	8 well Substrate Slides, Monkey Thyroid
1 x 0.5 ml	<b>CONTROL + ATA</b> *	Anti-Thyroid Positive Control. Contains human serum.
1 x 0.5 ml	<b>CONTROL -</b> *	Negative Control. Contains human serum.
1 x 5 ml	<b>IgG-CONJ FITC</b> *	Anti-human IgG FITC Conjugate. Protect from light.
1 x 5 ml	<b>IgG-CONJ FITC EB</b> **	Anti-human IgG FITC Conjugate containing Evan's Blue. Protect from light.
1 x 60 ml	<b>BUF</b> *	Buffered Diluent.
2 vials	<b>BUF WASH</b>	Phosphate Buffered Saline (PBS). Dissolve each vial to 1 liter.

- 1 x 5.0 ml **MOUNTING MEDIUM**\* Mounting Medium. Do not freeze.
- 1 x 1.0 ml **EVANS** Evan's Blue Counterstain.
- 1 x 12 **COVER SLD** Coverslips.

\* Contains < 0.1% NaN<sub>3</sub>

† Replaces conjugate without counterstain in Code numbers containing "EB"

**Material required but not provided**

- Fluorescence microscope
- Micropipette or Pasteur pipette
- Serological pipettes
- Staining dish (e.g. Coplin jar)
- Small test tubes (e.g. 13 x 75 mm) and test tube rack
- Distilled or deionized water
- 1 liter container
- Wash bottle
- Paper towels
- Incubation chamber

**WARNINGS AND PRECAUTIONS**

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>22</sup>.

WARNING - Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from IMMCO. Do not use beyond expiration date.

**SPECIMEN COLLECTION AND PREPARATION**

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

**PROCEDURE**

**Test Method**

**A. Screening**

1. Dilute each patient serum 1:10 with the Buffered Diluent provided (10 µl serum + 90 µl Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 50 µl of the Negative Control to well #1. Similarly apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 50 µl) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.

7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 50 µl) to each well.
9. Repeat steps **7 and 8** for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If desired, 2-3 drops of Evans blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.
12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

## B. Endpoint Determination (titration)

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Each test run should include the Positive and Negative Controls. Make serial two-fold dilutions starting at 1:10. The reciprocal of the highest dilution producing a positive reaction is the titer.

### Preparation of Serial Dilutions

Number four tubes 1 through 8. Add 0.9 ml of Sample Diluent to tube 1 and 0.2 ml to tubes 2 through 8. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

Tubes	1	2	3	4	5	6	7	8
Serum	0.1 ml							
	+							
Buffered Diluent	0.9 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
		↗	↗	↗	↗	↗	↗	↗
Transfer	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
Final dilution	1:10	1:20	1:40	1:80	1:160	1:320	1:640	1:1280 etc.

### QUALITY CONTROL

Both a positive and negative control serum should be included with each test run. The negative control should show no specific fluorescence, whereas the positive control should have 2+ or greater staining intensity .

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Gross contamination as a result of improper storage or handling. If signs of contamination such as turbidity are seen, discard and use another control.

- Problems with the optical system of the fluorescence microscope. These may include: improper alignment, use of the bulb beyond the expected performance life, etc.
- Allowing the slide to dry during the procedure.

### **INTERPRETATION OF RESULTS**

The results of the tests for thyroid antibodies should be read as negative (<10) or positive with titer. Read for specific staining of the colloid or the epithelial cells of the follicles. Various other tissue antibodies such as anti-nuclear antibodies (ANA), and anti-mitochondrial antibodies (AMA) may also be observed on thyroid tissue. Any sera giving nuclear staining reactions may be tested on Antinuclear Antibody (ANA) Test (HEp-2 Cells), Antinuclear Antibody (ANA) Test (Mouse Liver Sections). Any sera giving smooth muscle or mitochondrial staining reactions may be tested on Autoantibody Test System (Mouse Kidney/Stomach Sections).

### **LIMITATION OF THE PROCEDURE**

In some cases, sera positive for thyroid antibodies may either be very weak or negative at the initial screening dilution (prozone phenomenon). In such doubtful cases the sera should be screened at higher dilutions and, if positive, antibody titers determined.

In some cases the presence of two or more autoantibodies in a serum which are reactive with the same substrate may cause an interference in their detection by immunofluorescence. This interference may cause either a failure to detect anti-thyroid antibodies or a suppression of its titer if the interfering antibody has a higher titer.

### **CLINICAL SIGNIFICANCE**

Anti-thyroid antibodies are found in patients with autoimmune thyroid disorders and occasionally in healthy individuals. Table 1 at the end of this document summarizes the association of thyroid antibodies to the disease state<sup>3</sup>.

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### **REFERENCES•REFERENCIAS•LITERATUR•RIFERIMENTI**

1. Mooij P, Drexhage HA. Autoimmune thyroid disease. Clinics Lab Med. 13:683-697, 1993.
2. Champion BR, Cooke A, Reyner DC. Thyroid autoimmunity. Curr Opinion Immunol. 4:770-778, 1992.
3. Lindstedt G, Berg G, Jansson S, et al Clinical use of laboratory thyroid tests and investigations. J int Fed clin chem 6:136-141, 1994.
4. Pashke R, Vogg M, Swillens S, Usadel KH. Correlation of microsomal antibodies with the intensity of the intrathyroidal autoimmune process in Grave's disease. J Clin endocrinol Metab. 77:939-943, 1993.
5. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 3-40, 1987.
6. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1988 [HHS Pub. No. (CDC) 88-8395].

**Table 1: Significance of Thyroid Antibodies**

Ab Specificity	Disease Association	Indications for use
<b>Thyroglobulin</b>	1) thyroid autoimmune diseases, 2) increased concentrations after thyroidectomy for undifferentiated carcinoma indicates metastasis 3) incorrect thyroglobulin levels.	1) complement to thyroglobulin levels 2) follow up of differentiated thyroid carcinoma
<b>Microsomal (TPO)</b>	Thyroid autoimmune disease	1) TSH rise of unknown cause 2) Goitre of unknown etiology 3) Evaluation of suspected autoimmune polyglandular endocrine disorders. 4) Familial Studies. 5) Assessment of risk of thyroid dysfunction during drug treatment affecting the thyroid, and/or the immune system. 6) Screening for risk of postpartum thyroiditis 7) Hyperthyroidism

# NOTES

# NOTES

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**

**60 Pineview Drive**

**Buffalo, NY 14228-2120**

**Telephone: (716) 691-0091**

**Fax: (716) 691-0466**

**Toll Free USA/Canada: 1-800-537-TEST**

**E-Mail: [info@immcodiagnostics.com](mailto:info@immcodiagnostics.com)**

*or your local product distributor*

EU Authorized Representative/Autorisierter Repräsentant/Rappresentante  
Autorizzato/Representante Autorizado/Représentant Autorisé

EMERGO Group, Inc.

Molenstraat 15, 2513 BH, The Hague,

The Netherlands

Tel (+31) 345 8570, Fax (+31) 346 7299

[www.emergogroup.com](http://www.emergogroup.com)